HYALOFAST Čeština

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2016/07

HYALOFAST English

Support for human mesenchymal stem cells

Sterilized by radiation.

STERILE | R | Don't use if package is opened or damaged

PRODUCT INFORMATION AND INDICATIONS

Mesenchymal stem cells (MSCs) residing in bone marrow are the progenitors for chondrocytes, osteoblasts and several other types of cells.

HYALOFAST is a biologically biodegradable support for the entrapment of mesenchymal stem cells for the repair of chondral or osteochondral lesions.

HYALOFAST is initially filling the defect and temporarily substituting the chondral and osteochondral tissue in the lesions until it is absorbed and replaced by the skeletal tissue. HYALOFAST may act as a support for bone marrow aspirate or as a chondroprotective coverage which favours in situ residence of mesenchymal stem cells after their mobilization due to microfracture or perforation procedures.

HYALOFAST is a white non-woven pad entirely composed of HYAFF®, a semi-synthetic derivative of hyaluronic acid, a naturally occurring extracellular matrix molecule and a major component of human cartilage.

The main characteristics of HYALOFAST are the following:

- complete biocompatibility of the support: the polymer used is an ester derivative of hyaluronic acid (one of the main components of the extracellular matrix). After spontaneous degradation, the polymer releases the parent molecule enriching the transplantation site with hyaluronic acid.
- three-dimensional structure: the non-woven mesh favours cell adhesion, therefore MSCs are able to arrange themselves inside the artificial three-dimensional matrix.
- ease of handling and application
- the implant is soft and conformable and readily adheres to the site of application

WARNINGS AND PRECAUTIONS

- HYALOFAST is for exclusive use in surgical procedures.
- HYALOFAST should be used immediately after opening of the pouch.
- HYALOFAST is sterilized by gamma-radiation and is for single use only. Portions of unutilized product must be discarded.

- There is a risk of non sterile product and patient infection, in case any unutilized portions of product were to be re-used after first application.
- Sterility is guaranteed as long as the package is closed and undamaged.
- Do not use in case of package damage and report to local distributor.
- Store in a cool and dry place (Ts 40°C).
- After use, dispose according to applicable national practice.
- HYALOFAST does not provide a structural support after implantation, therefore the use of a proper load bearing device is recommended, if needed.

KONTRAINDICATIONS

The use of HYALOFAST is contraindicated in patients with hypersensitivity to the product components.

INSTRUCTIONS FOR USE

Prior to use, the product's secondary packaging is removed and HYALOFAST is transported in its primary packaging to the sterile field. Both sides of the scaffold are equivalent for application.

HYALOFAST can be:

- pre-loaded with bone marrow aspirate prior to implantation;
- applied directly onto the debrided areas after microfracture or perforation procedures, with the aim of favouring in situ residence of cells mobilized by such procedures.

For both approaches, the recommended surgical procedure for application of the product is the following. The implantation of HYALOFAST can be performed in miniarthrotomy or arthroscopy, at the physician discretion. After joint inspection, the chondral or osteochondral lesion is carefully debrided until a stable shoulder surrounds the defect. Defect size can be assessed by a sterile template.

HYALOFAST is cut with sterile scissors to fit the defect area and placed onto the lesion with forceps.

When applied in arthroscopy, after creating a dry joint surface HYALOFAST is positioned into the lesion with the help of a flattened probe.

If necessary, HYALOFAST may be secured to the defect by standard fixation means such as fibrin glue, sutures or pins. In presence of very deep defects, the overlapping of several pads could be used as long as the total thickness of the implant is not greater than the depth of the defect.

Before and after tourniquet release, the joint is mobilized to ensure adherence and stability prior to standard closure of the surgical incision.

At the surgeon's discretion an articular drainage and a compressive bandage may be applied for 24 hours post-implantation. Subsequent immobilization is recommended for 24 hours post-surgery, after which a postoperative rehabilitation program can be initiated.

PRESENTATIONS

HYALOFAST is available in different presentations:

- pouch containing 1 pad 2x2 cm
- pouch containing 1 pad 5x5 cm

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Anika Therapeutics S.r.l.