

HemiCAP

PF W



HemiCAP
Patello-Femoral Arthroplasty

Sistema Hemicap PF W (Referencias)

System Catalog

Instrumentation System

7000-2300	Instrument Kit, Patello-Femoral W
7000-2302	Instrument Kit, 25/30mm Patella
PX07-1205	2.5mm Guide Pin, W (5 Pk) (non - sterile)
PX00-0200	2.5mm Guide Pin Kit, W (sterile)

Articular Component, Patella

P255-1050	25mm Anatomic, 7.0mm thick
P306-0070	30mm Button, 7.0mm thick
P306-0090	30mm Dome, 9.0mm thick

Taper Post

PX11-0218	Taper Post, 11mm x 21.5mm
-----------	---------------------------

Fixation Stud

PX75-0173	Fixation Stud, 7.5mm x 18.5mm
-----------	-------------------------------

Articular Component, Trochlear

	⊕ S/I	⊖ M/L
PX02-0704	7.0mm x 4.0mm Offset	
PX02-0705	7.0mm x 5.0mm Offset	
PX02-0854	8.5mm x 4.0mm Offset	
PX02-0855	8.5mm x 5.0mm Offset	
PX02-1004	10.0mm x 4.0mm Offset	
PX02-1005	10.0mm x 5.0mm Offset	
PX02-1154	11.5mm x 4.0mm Offset	
PX02-1155	11.5mm x 5.0mm Offset	

Diane L. Dahm, MD. Patellofemoral Arthroplasty Versus Total Knee Arthroplasty in Patients with Isolated Patellofemoral Osteoarthritis. The American Journal of Orthopedics, October 2010

Materials

Articular Component:

Cobalt-Chromium Alloy (Co-Cr-Mo)

Surface Coating:

Titanium (CPTi)

Taper Post:

Titanium Alloy (Ti-6Al-4V)

Tibial Component:

Ultra-High-Molecular Weight Polyethylene (UHMWPE)

Fabricado por:


arthrosurface®

HemiCAP

PF HemiCAP®
KAHUNA



HemiCAP *Kahuna*
Patello-Femoral Arthroplasty

Sistema Hemicap PF KAHUNA
(Referencias)

System Catalog

Articular Component, Trochlear

PX11-0218	PF W Poste de fijación 11 mm
PX75-0173	PF Poste de fijación 7.5 mm
PWL2-0705	7.0 x 5.0 mm, Izq.
PWL2-0855	8.5 x 5.0 mm, Izq.
PWL2-1005	10.0 x 5.0 mm, Izq.
PWR2-0705	7.0 x 5.0 mm, drch.
PWR2-0855	8.5 x 5.0 mm, drch.
PWR2-1005	10.0 x 5.0 mm, drch.
PW09-1500	Fresa (fungible)

Materials

Articular Component:

Surface Coating:

Taper Post:

Tibial Component:

Cobalt-Chromium Alloy (Co-Cr-Mo)

Titanium (CPTi)

Titanium Alloy (Ti-6Al-4V)



Ultra-High-Molecular Weight Polyethylene (UHMWPE)

Fabricado por:


arthrosurface®

INSTRUMENTACIÓN


Hemicap® PF Classic




INSTRUMENTACIÓN

Hemicap® PF W

BANDEJA SUPERIOR



BANDEJA INFERIOR



INSTRUMENTAL PATELA

